Size:165x120mm

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BOUILLOTTE
MICRO-ONDES

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EN
MANUAL
As a hot compress:
Remove the cushion from its cover before reheating it in the microwave.
Check that the turntable is clean. Place the cushion on the tray, distributing its weight evenly. Add a half-filled cup of water to the microwave. Heat for up to 1 minute, checking the temperature and rotating the pad at 30 second intervals to prevent it from burning. Refer to the table to know the heating time according to the power of your microwave. Do not heat unattended. Take out of the microwave, let the cushion rest for 20 to 30 seconds. Always cover the cushion with its cover before applying it to the body.

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HEAT AGAIN ONLY WHEN ITEM HAS RETURNED TO A ROOM TEMPERATURE.

As a cold compress:
Place the product in a plastic bag and then in the freezer until the desired temperature is reached.

UPKEEP
The cover can be hand or machine washed cold, protected in a canvas bag. Wipe the product with a damp cloth and allow to dry before use.

CAUTION
Carefully examine the cushion for damage before using it. Always remove the cover before placing the cushion in the microwave. Always cover the cushion with its cover before using it to avoid direct contact of heat with the skin. Do not overheat or exceed the one minute heating limit. Do not ingest content. If product is too hot, set aside and allow to cool before use. Never apply a hot or cold compress to irritated or damaged skin. If you feel any discomfort or pain discontinue use immediately.
Not suitable for children under 3 years old.

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